



PROCESSES

Sterilization
Washing
Chemical
bio-decontamination
Contamination control
Dense-phase fluids

SPECIAL APPLICATIONS

Rotating sterilizers
Pharmaceutical
closures
processing

COMPONENTS

Pressure vessels
Doors
Piping
Process controller
Wireless sensors

ACCESSORIES

Un/loading conveyors
Customized racks

Multi-Purpose Sterilizers

Loading the closures



1. Closures are manually loaded into the rotating drum.

2. The drum is loaded into the autoclave chamber.

Docking the bagging cart



3. Closures are washed, siliconized, sterilized and dried in one processes. At the end the load is discharged without any contact with the operator. Various options for discharging or transfer are available to meet the most different integration requirements.

Pharmaceutical Elastomeric Closures

CONTACT US:

Italy

Fedegari Autoclavi SpA
SS 235 km 8
27100 Albuzzano (PV)
T +39 0382 434111
F +39 0382 434150
E fedegari@fedegari.com

Germany

Fedegari GmbH
Lehrer-Götz-Weg, 11
81825 Munich
T +49 89 42044920
F +49 89 42044918
E info@de.fedegari.com

Singapore

Fedegari Asia Pte. Ltd.
8, Boon Lay Way
#03-16 - Tradehub 21
609964 Singapore
T +65 6 3164761
F +65 6 8967365
E info@fedegariasia.com


Switzerland

Fedegari (SUISSE) SA
Via alla Gerra, 11
6930 Bedano
T +41 91 9352090
F +41 91 9352099
E info@fedegarisuisse.ch

USA

Fedegari Technologies, Inc.
1228 Bethlehem Pike
18960 Sellersville (PA)
T +1 215 4530400
F +1 215 4530406
E info@fedegariusa.com

Discharging the closures in a continuous PE liner



4. Slide and rotate the adapter housing.



5. Install the liner on the adapter and pull it through the welding brackets.



6. After welding, pull the liner all the way down until it reaches the scale.



7. The closures are discharged inside the liner. When the desired volume is reached the interface stops the drum inside the chamber. Weld the liner and repeat the operation until the batch has been bagged.